



CDK14 Polyclonal Antibody

Cat No: HR1AP11321

For research use only

Overview

Product Name	CDK14 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	51kDa
GeneID?Human?	CDK14 KIAA0834 PFTK1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	cyclin dependent kinase 14(CDK14) Homo sapiens PFTK1 is a member of the CDC2 (MIM 116940)-related protein kinase family (Yang and Chen, 2001 [PubMed 11313143]).[supplied by OMIM, Mar 2008],