

RL1D1 Polyclonal Antibody

Cat No: HR1AP11291

For research use only

Overview

Product Name	RL1D1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	53kDa
GeneID?Human?	RSL1D1 CATX11 CSIG PBK1 L12
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the ribosomal protein L1P family. Highly divergent.,tissue specificity:Placenta.,