



## Serinc1 Polyclonal Antibody

Cat No: HR1AP4540

For research use only

### Overview

Product Name	Serinc1 Polyclonal Antibody	
Source	Rabbit	
Applications	WB,IHC-p,IF,ELISA	
Species Reactivity	Human,Mouse,Rat	
Recommended Dilutions		
Immunogen		
Species	Rabbit	
Storage	-20°C/1 year	
Isotype		
Clonality		
Concentration	1 mg/ml	
Observed band	48kDa	
GenelID?Human?	SERINC1	
Human Swiss-Prot No.		
Cellular localization		
Alternative Names	SERINC1; KIAA1253; TDE1L; TDE2; Serine incorporator 1; Tumor differentially expressed protein 1-like; Tumor differentially expressed protein 2	
Background	function:Enhances the incorporation of serine into phosphatidylserine and sphingolipids.,similarity:Belongs to the TDE1 family.,subunit:Interacts with SPTLC1.,	