

SGK1 (4D12) Mouse mAb

Cat No: HR1AM2375

For research use only

Overview

Product Name	SGK1 (4D12) Mouse mAb
Source	Mouse
Applications	IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	IHC 1:100-200
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	45-60kDa
GeneID?Human?	6446
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SGK variant i3, serum/glucocorticoid regulated kinase 1
Background	Serum and glucocorticoid-inducible kinase (SGK) is a serine/threonine kinase closely related to Akt. SGK is rapidly induced in response to a variety of stimuli, including serum, glucocorticoid, follicle stimulating hormone, osmotic shock, and mineralocorticoids.