



UNC5A Polyclonal Antibody

Cat No: HR1AP11001

For research use only

Overview

Product Name	UNC5A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	92kDa
GeneID?Human?	UNC5A KIAA1976 UNC5H1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	unc-5 netrin receptor A(UNC5A) Homo sapiens UNC5A belongs to a family of netrin-1 (MIM 601614) receptors thought to mediate the chemorepulsive effect of netrin-1 on specific axons. For more information on UNC5 proteins, see UNC5C (MIM 603610).[supplied by OMIM, Apr 2004],