

## IL-8 (13F8) Mouse mAb Cat No: HR1AM2371

For research use only

## Overview

Product Name	IL-8 (13F8) Mouse mAb
Source	Mouse
Applications	IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	IHC 1:100-200
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	11kDakDa
GeneID?Human?	3576
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	9E3,CXCL8,CEF-4,GCP-1,Granulocyte chemotactic protein 1
Background	IL-8 is a chemotactic factor that attracts neutrophils, basophils, and T-cells, but not monocytes. It is also involved in neutrophil activation. It is released from several cell types in response to an inflammatory stimulus.