

TF2LX Polyclonal Antibody

Cat No: HR1AP10854

For research use only

Overview

Product Name	TF2LX Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	26kDa
GeneID?Human?	TGIF2LX TGIFLX
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	TGFB induced factor homeobox 2 like, X-linked(TGIF2LX) Homo sapiens This gene encodes a member of the TALE/TGIF homeobox family of transcription factors. Testis-specific expression suggests that this gene may play a role in spermatogenesis. A homolog of this gene lies within the male specific region of chromosome Y, in a block of sequence that is thought to be the result of a large X-to-Y transposition. [provided by RefSeq, Jul 2008],