

BNIP-3 Polyclonal Antibody

Cat No: HR1AP8363

For research use only

Overview

Product Name	BNIP-3 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	BNIP3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	BNIP3; NIP3; BCL2/adenovirus E1B 19 kDa protein-interacting protein 3
Background	BCL2 interacting protein 3(BNIP3) Homo sapiens This gene is encodes a mitochondrial protein that contains a BH3 domain and acts as a pro-apoptotic factor. The encoded protein interacts with anti-apoptotic proteins, including the E1B 19 kDa protein and Bcl2. This gene is silenced in tumors by DNA methylation. [provided by RefSeq, Dec 2014],