

STAT1 (12E8) Mouse mAb

Cat No: HR1AM2370

For research use only

Overview

Product Name	STAT1 (12E8) Mouse mAb
Source	Mouse
Applications	IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	IHC 1:100-200
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	84,91kDakDa
GenelD?Human?	6772
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CANDF7 antibody, ISGF3 antibody, DKFZp686B04100 antibody
Background	STAT1 is a member of the Signal Transducers and Activators of Transcription family of transcription factors. STAT1 is involved in upregulating genes due to a signal by either type II, type III, or type III interferons.