

MARH7 Polyclonal Antibody

Cat No: HR1AP10739

For research use only

Overview

Product Name	MARH7 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	77kDa
GeneID?Human?	MARCH7 AXOT RNF177
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	membrane associated ring-CH-type finger 7(MARCH7) Homo sapiens MARCH7 is a member of the MARCH family of membrane-bound E3 ubiquitin ligases (EC 6.3.2.19). MARCH proteins add ubiquitin (see MIM 191339) to target lysines in substrate proteins, thereby signaling their vesicular transport between membrane compartments (Bartee et al., 2004 [PubMed 14722266]).[supplied by OMIM, Mar 2010],