JAK1 (7G6) Mouse mAb

Cat No: HR1AM2368

For research use only

Overview

Product Name	JAK1 (7G6) Mouse mAb
Source	Mouse
Applications	IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	IHC 1:100-200
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	130kDakDa
GeneID?Human?	3716
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Janus kinase 1,JTK3,Tyrosine protein kinase JAK1,JAK1A,JAK1B
Background	Members of the Janus family of tyrosine kinases (Jak1, Jak2, Jak3 and Tyk2) are activated by ligands binding to a number of associated cytokine receptors.