

ARP19 Polyclonal Antibody

Cat No: HR1AP10600

For research use only

Overview

Product Name	ARP19 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	12kDa
GeneID?Human?	ARPP19
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	cAMP regulated phosphoprotein 19(ARPP19) Homo sapiens The 19-kD cAMP-regulated phosphoprotein plays a role in regulating mitosis by inhibiting protein phosphatase-2A (PP2A; see MIM 176915) (summary by Gharbi-Ayachi et al., 2010 [PubMed 21164014]).[supplied by OMIM, Feb 2011],