

CAB39 Polyclonal Antibody

Cat No: HR1AP10582

For research use only

Overview

Product Name	CAB39 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	CAB39 MO25 CGI-66
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Together with the STE20-related adaptor-alpha (STRAD alpha) pseudo kinase, forms a regulatory complex capable of stimulating the activity of STK11.,similarity:Belongs to the Mo25 family.,