

B2L14 Polyclonal Antibody

Cat No: HR1AP10510

For research use only

Overview

Product Name	B2L14 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	35kDa
GeneID?Human?	BCL2L14 BCLG
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>BCL2 like 14(BCL2L14) Homo sapiens The protein encoded by this gene belongs to the BCL2 protein family. BCL2 family members form hetero- or homodimers and act as anti- or pro-apoptotic regulators that are involved in a wide variety of cellular activities. Overexpression of this gene has been shown to induce apoptosis in cells. Three alternatively spliced transcript variants encoding two distinct isoforms have been reported for this gene. [provided by RefSeq, May 2009],</p>