

B3A2 Polyclonal Antibody

Cat No: HR1AP10504

For research use only

Overview

Product Name	B3A2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	136kDa
GeneID?Human?	SLC4A2 AE2 EPB3L1 HKB3 MPB3L
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	alternative products:Additional isoforms seem to exist,function:Plasma membrane anion exchange protein of wide distribution.,similarity:Belongs to the anion exchanger (TC 2.A.31) family.,