

# IFN14 Polyclonal Antibody

**Cat No: HR1AP10273**

For research use only

## Overview

|                       |                                                                                                                                                                                                                                         |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | IFN14 Polyclonal Antibody                                                                                                                                                                                                               |
| Source                | Rabbit                                                                                                                                                                                                                                  |
| Applications          | WB,ELISA                                                                                                                                                                                                                                |
| Species Reactivity    | Human                                                                                                                                                                                                                                   |
| Recommended Dilutions |                                                                                                                                                                                                                                         |
| Immunogen             |                                                                                                                                                                                                                                         |
| Species               | Rabbit                                                                                                                                                                                                                                  |
| Storage               | -20°C/1 year                                                                                                                                                                                                                            |
| Isotype               |                                                                                                                                                                                                                                         |
| Clonality             |                                                                                                                                                                                                                                         |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                 |
| Observed band         | 20kDa                                                                                                                                                                                                                                   |
| GeneID?Human?         | IFNA14                                                                                                                                                                                                                                  |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                         |
| Cellular localization |                                                                                                                                                                                                                                         |
| Alternative Names     |                                                                                                                                                                                                                                         |
| Background            | function:Produced by macrophages, IFN-alpha have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase.,similarity:Belongs to the alpha/beta interferon family., |