



## Lamin B1 Polyclonal Antibody

Cat No: HR1AP5021

For research use only

### Overview

Product Name	Lamin B1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	70kDa
GenID?Human?	LMNB1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	LMNB1; LMN2; LMNB; Lamin-B1
Background	lamin B1(LMNB1) <i>Homo sapiens</i> This gene encodes one of the two B-type lamin proteins and is a component of the nuclear lamina. A duplication of this gene is associated with autosomal dominant adult-onset leukodystrophy (ADLD). Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2015],