

## SRBP2 Polyclonal Antibody Cat No: HR1AP10094

For research use only

## Overview

Product Name	SRBP2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	125kDa
GeneID?Human?	SREBF2 BHLHD2 SREBP2
Human Swiss-Prot No.	
Cellular	
localization	
Alternative Names	
Background	sterol regulatory element binding transcription factor 2(SREBF2) Homo sapiens This gene encodes a member of the a ubiquitously expressed transcription factor that controls cholesterol homeostasis by regulating transcription of sterol-regulated genes. The encoded protein contains a basic helix-loop-helix-leucine zipper (bHLH-Zip) domain and binds the sterol regulatory element 1 motif. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013],