

I κ B β ? Mouse mAb

Cat No: HR1AM2337

For research use only

Overview

Product Name	I κ B β ? Mouse mAb
Source	Mouse
Applications	IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	IHC 1:100-200
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	48kDa
GeneID?Human?	4793
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	I κ B β , I κ B β , NF κ B β , TRIP9
Background	The NF- κ B/Rel transcription factors are present in the cytosol in an inactive state complexed with the inhibitory I κ B proteins.