

HSPB8/HSP22 Mouse mAb

Cat No: HR1AM2341

For research use only

Overview

Product Name	HSPB8/HSP22 Mouse mAb
Source	Mouse
Applications	WB
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:500-1,000
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	22kDakDa
GeneID?Human?	26353
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Alpha crystallin C chain antibody, Heat shock 22KDa protein 8 antibody, HMN 2 antibody
Background	Heat shock protein beta-8 is a protein that in humans is encoded by the HSPB8 gene. HSPB8 (HSP22) is a member of the small heat shock protein superfamily and the human protein is most closely related to HSP27.