



AMPK α 1 Mouse mAb

Cat No: HR1AM2336

For research use only

Overview

Product Name	AMPK α 1 Mouse mAb
Source	Mouse
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:100-200
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	63kDa
GeneID?Human?	5562
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	AAPK1, α 1 AMP activated protein kinase catalytic subunit, PRKAA1, SNF1A, MGC33776
Background	AMP-activated protein kinase (AMPK) is highly conserved from yeast to plants and animals and plays a key role in the regulation of energy homeostasis. AMPK is a heterotrimeric complex composed of a catalytic α subunit and regulatory β and γ subunits, each of which is encoded by two or three distinct genes (β 1, β 2; γ 1, γ 2; γ 1, γ 2, γ 3).