

MOT4 Polyclonal Antibody

Cat No: HR1AP10041

For research use only

Overview

Product Name	MOT4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	50kDa
GeneID?Human?	SLC16A3 MCT4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Monocarboxylate transporter 4 (MCT 4) (Solute carrier family 16 member 3)
Background	<p>solute carrier family 16 member 3(SLC16A3) Homo sapiens Lactic acid and pyruvate transport across plasma membranes is catalyzed by members of the proton-linked monocarboxylate transporter (MCT) family, which has been designated solute carrier family-16. Each MCT appears to have slightly different substrate and inhibitor specificities and transport kinetics, which are related to the metabolic requirements of the tissues in which it is found. The MCTs, which include MCT1 (SLC16A1; MIM 600682) and MCT2 (SLC16A7; MIM 603654), are characterized by 12 predicted transmembrane domains (Price et al., 1998 [PubMed 9425115]).[supplied by OMIM, Mar 2008].</p>