

TIF-IA Polyclonal Antibody

Cat No: HR1AP10023

For research use only

Overview

| | |
|-----------------------|---|
| Product Name | TIF-IA Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human,Mouse,Rat |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 70kDa |
| GeneID?Human? | RRN3 TIFIA |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | RNA polymerase I-specific transcription initiation factor RRN3 (Transcription initiation factor IA) (TIF-IA) |
| Background | function:Required for efficient transcription initiation by RNA polymerase I.,similarity:Belongs to the RRN3 family., |