

ATG13 Rabbit Polyclonal Antibody

Cat No: HR1AP9685

For research use only

Overview

Product Name	ATG13 Rabbit Polyclonal Antibody
Source	Rabbit
Applications	IHC-p
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	57kDa
GeneID?Human?	ATG13 KIAA0652
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Autophagy-related protein 13
Background	autophagy related 13(ATG13) Homo sapiens The protein encoded by this gene is an autophagy factor and a target of the TOR kinase signaling pathway. The encoded protein is essential for autophagosome formation and mitophagy. [provided by RefSeq, Oct 2016],