

## TBP/TATA Binding Protein Mouse mAb

Cat No: HR1AM2331

For research use only

## **Overview**

Product Name	TBP/TATA Binding Protein Mouse mAb
1 Toddet Ivallie	TEL / TATA Billiding Flotelin Modes HIAD
Source	Mouse
Applications	IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	IHC 1:100-200
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	35-45kDakDa
GeneID?Human?	6908
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GFT2D antibody, MGC117320 antibody, SCA17 antibody, TATA binding protein antibody, TF2D antibody
Background	The TATA-binding protein (TBP) is a general transcription factor that binds specifically to a DNA sequence called the TATA box. This DNA sequence is found about 30 base pairs upstream of the transcription start site in some eukaryotic gene promoters.