

Dab1 (phospho Tyr232) Polyclonal Antibody

Cat No: HR1AP8963

For research use only

Overview

Product Name	Dab1 (phospho Tyr232) Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	60kDa
GeneID?Human?	DAB1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	DAB1; Disabled homolog 1
Background	DAB1, reelin adaptor protein(DAB1) Homo sapiens The laminar organization of multiple neuronal types in the cerebral cortex is required for normal cognitive function. In mice, the disabled-1 gene plays a central role in brain development, directing the migration of cortical neurons past previously formed neurons to reach their proper layer. This gene is similar to disabled-1, and the protein encoded by this gene is thought to be a signal transducer that interacts with protein kinase pathways to regulate neuronal positioning in the developing brain. Alternatively spliced transcript variants of this gene have been reported, but their full length nature has not been determined. [provided by RefSeq, Jul 2008],