

Desmin Mouse mAb

Cat No: HR1AM2329

For research use only

Overview

Product Name	Desmin Mouse mAb
Source	Mouse
Applications	IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	IHC 1:100-200
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	53kDakDa
GeneID?Human?	1674
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CMD11, CSM1,CSM2, DES, FLJ12025
Background	Desmin is a myogenic marker expressed in early development that forms a network of filaments that extends across the myofibril and surrounds Z discs. The desmin cytoskeleton provides a connection among myofibrils, organelles and the cytoskeleton.