



## Glut1 Polyclonal Antibody

Cat No: HR1AP3549

For research use only

### Overview

Product Name	Glut1 Polyclonal Antibody
Source	Rabbit
Applications	IF,WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	55kDa
GenelD?Human?	SLC2A1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SLC2A1; GLUT1; Solute carrier family 2; facilitated glucose transporter member 1; Glucose transporter type 1, erythrocyte/brain; GLUT-1; HepG2 glucose transporter
Background	solute carrier family 2 member 1(SLC2A1) <i>Homo sapiens</i> This gene encodes a major glucose transporter in the mammalian blood-brain barrier. The encoded protein is found primarily in the cell membrane and on the cell surface, where it can also function as a receptor for human T-cell leukemia virus (HTLV) I and II. Mutations in this gene have been found in a family with paroxysmal exertion-induced dyskinesia. [provided by RefSeq, Apr 2013],