



Netrin-1 Polyclonal Antibody

Cat No: HR1AP4032

For research use only

Overview

Product Name	Netrin-1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	75kDa
GeneID?Human?	NTN1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NTN1; NTN1L; Netrin-1; Epididymis tissue protein Li 131P
Background	netrin 1(NTN1) Homo sapiens Netrin is included in a family of laminin-related secreted proteins. The function of this gene has not yet been defined; however, netrin is thought to be involved in axon guidance and cell migration during development. Mutations and loss of expression of netrin suggest that variation in netrin may be involved in cancer development. [provided by RefSeq, Jul 2008],