



TRAPPC5 Polyclonal Antibody

Cat No: HR1AP5829

For research use only

Overview

Product Name	TRAPPC5 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	TRAPPC5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TRAPPC5; Trafficking protein particle complex subunit 5
Background	function:May play a role in vesicular transport from endoplasmic reticulum to Golgi.,similarity:Belongs to the TRAPP small subunits family. BET3 subfamily.,subunit:Part of the multisubunit TRAPP (transport protein particle) complex.,