

# Cyclophilin B Mouse mAb

Cat No: HR1AM2330

For research use only

## Overview

Product Name	Cyclophilin B Mouse mAb
Source	Mouse
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:100-200
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	21kDa
GeneID?Human?	5479
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Cphn2, MGC2224, Rotamase, S cyclophilin, CYPB, CyP-S1, Rotamase B
Background	Cyclophilin B (CyPB) is an ER-localized chaperone protein belonging to the family of peptidyl-prolyl cis-trans isomerases.