




SMVT Polyclonal Antibody

Cat No: HR1AP9012

For research use only

Overview

Product Name	SMVT Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	68kDa
GenelD?Human?	SLC5A6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SLC5A6; SMVT; Sodium-dependent multivitamin transporter; Na(+) -dependent multivitamin transporter; Solute carrier family 5 member 6
Background	function: Transports pantothenate, biotin and lipoate in the presence of sodium., similarity: Belongs to the sodium:solute symporter (SSF) (TC 2.A.21) family.,



bioelsa