




SFT2D3 Polyclonal Antibody

Cat No: HR1AP8892

For research use only

Overview

Product Name	SFT2D3 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GenID?Human?	SFT2D3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SFT2D3; Vesicle transport protein SFT2C; SFT2 domain-containing protein 3
Background	function:May be involved in fusion of retrograde transport vesicles derived from an endocytic compartment with the Golgi complex.,similarity:Belongs to the SFT2 family.,



bioelsa