




SFT2D2 Polyclonal Antibody

Cat No: HR1AP7108

For research use only

Overview

Product Name	SFT2D2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	18kDa
GenelID?Human?	SFT2D2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SFT2D2; Vesicle transport protein SFT2B; SFT2 domain-containing protein 2
Background	function:May be involved in fusion of retrograde transport vesicles derived from an endocytic compartment with the Golgi complex.,similarity:Belongs to the SFT2 family.,



bioelsa