

GPR174 Polyclonal Antibody

Cat No: HR1AP8880

For research use only

Overview

Product Name	GPR174 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	34kDa
GeneID?Human?	GPR174
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPR174; FKSG79; Probable G-protein coupled receptor 174
Background	G protein-coupled receptor 174(GPR174) Homo sapiens This gene encodes a protein belonging to the G protein-coupled receptor superfamily. These proteins are characterized by the presence of seven alpha-helical transmembrane domains, and they activate or interact with various endogenous or exogenous ligands, including neurotransmitters, hormones, and odorant and taste substances. This family member is classified as an orphan receptor because the cognate ligand has not been identified. [provided by RefSeq, Sep 2011],