

## GPR142 Polyclonal Antibody

Cat No: HR1AP6972

For research use only

### Overview

Product Name	GPR142 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	51kDa
GeneID?Human?	GPR142
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPR142; PGR2; Probable G-protein coupled receptor 142; G-protein coupled receptor PGR2
Background	G protein-coupled receptor 142(GPR142) Homo sapiens GPR142 is a member of the rhodopsin family of G protein-coupled receptors (GPRs) (Fredriksson et al., 2003 [PubMed 14623098]).[supplied by OMIM, Mar 2008],