



Myoglobin Polyclonal Antibody

Cat No: HR1AP7309

For research use only

Overview

Product Name	Myoglobin Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	17kDa
GenelD?Human?	MB
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MB; Myoglobin
Background	myoglobin(MB) <i>Homo sapiens</i> This gene encodes a member of the globin superfamily and is expressed in skeletal and cardiac muscles. The encoded protein is a haemoprotein contributing to intracellular oxygen storage and transcellular facilitated diffusion of oxygen. At least three alternatively spliced transcript variants encoding the same protein have been reported. [provided by RefSeq, Jul 2008],