



# Frizzled-5/8 Polyclonal Antibody

Cat No: HR1AP8672

For research use only

## Overview

|                       |                                                                                                                                                                                                                                                                                |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | Frizzled-5/8 Polyclonal Antibody                                                                                                                                                                                                                                               |
| Source                | Rabbit                                                                                                                                                                                                                                                                         |
| Applications          | IHC-p,IF,ELISA                                                                                                                                                                                                                                                                 |
| Species Reactivity    | Human,Mouse,Rat                                                                                                                                                                                                                                                                |
| Recommended Dilutions |                                                                                                                                                                                                                                                                                |
| Immunogen             |                                                                                                                                                                                                                                                                                |
| Species               | Rabbit                                                                                                                                                                                                                                                                         |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                                                   |
| Isotype               |                                                                                                                                                                                                                                                                                |
| Clonality             |                                                                                                                                                                                                                                                                                |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                                        |
| Observed band         | kDa                                                                                                                                                                                                                                                                            |
| GeneID?Human?         | FZD5/FZD8                                                                                                                                                                                                                                                                      |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                                                |
| Cellular localization |                                                                                                                                                                                                                                                                                |
| Alternative Names     | FZD5; C2orf31; Frizzled-5; Fz-5; hFz5; FzE5; FZD8; Frizzled-8; Fz-8; hFz8                                                                                                                                                                                                      |
| Background            | frizzled class receptor 5(FZD5) Homo sapiens Members of the ‘frizzled’ gene family encode 7-transmembrane domain proteins that are receptors for Wnt signaling proteins. The FZD5 protein is believed to be the receptor for the Wnt5A ligand. [provided by RefSeq, Jul 2008], |