



GOLT1A Polyclonal Antibody

Cat No: HR1AP5841

For research use only

Overview

| | |
|-----------------------|---|
| Product Name | GOLT1A Polyclonal Antibody |
| Source | Rabbit |
| Applications | IHC-p,ELISA |
| Species Reactivity | Human,Mouse |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | kDa |
| GenID?Human? | GOLT1A |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | GOLT1A; GOT1B; Vesicle transport protein GOT1A; Golgi transport 1 homolog A; hGOT1b |
| Background | function:May be involved in fusion of ER-derived transport vesicles with the Golgi complex.,similarity:Belongs to the GOT1 family., |