



BRS-3 Polyclonal Antibody

Cat No: HR1AP8424

For research use only

Overview

Product Name	BRS-3 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	47kDa
GeneID?Human?	BRS3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	BRS3; Bombesin receptor subtype-3; BRS-3
Background	bombesin receptor subtype 3(BRS3) Homo sapiens The protein encoded by this gene is a G protein-coupled membrane receptor that binds bombesin-like peptides. This binding results in activation of a phosphatidylinositol-calcium second messenger system, with physiological effects including regulation of metabolic rate, glucose metabolism, and hypertension. [provided by RefSeq, Sep 2011],