



ERGIC-3 Polyclonal Antibody

Cat No: HR1AP7659

For research use only

Overview

| | |
|-----------------------|---|
| Product Name | ERGIC-3 Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,IHC-p,IF,ELISA |
| Species Reactivity | Human,Mouse |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 43kDa |
| GenID?Human? | ERGIC3 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | ERGIC3; C20orf47; ERV46; SDBCAG84; CGI-54; Endoplasmic reticulum-Golgi intermediate compartment protein 3; Serologically defined breast cancer antigen NY-BR-84 |
| Background | function:Possible role in transport between endoplasmic reticulum and Golgi.,similarity:Belongs to the ERGIC family.,subcellular location:Cycles between the endoplasmic reticulum and the Golgi.,subunit:Interacts with ERGIC1/ERGIC32., |