

# ApoL4 Polyclonal Antibody

Cat No: HR1AP8737

For research use only

## Overview

Product Name	ApoL4 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	34kDa
GeneID?Human?	APOL4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	APOL4; Apolipoprotein L4; Apolipoprotein L-IV; ApoL-IV
Background	apolipoprotein L4(APOL4) Homo sapiens The protein encoded by this gene is a member of the apolipoprotein L family and may play a role in lipid exchange and transport throughout the body, as well as in reverse cholesterol transport from peripheral cells to the liver. Two transcript variants encoding two different isoforms have been found for this gene. Only one of the isoforms appears to be a secreted protein. [provided by RefSeq, Jul 2008].