



## Arfaptin 1 Polyclonal Antibody

Cat No: HR1AP6657

For research use only

### Overview

|                       |   |
|-----------------------|---|
| Product Name          | Arfaptin 1 Polyclonal Antibody  |
| Source                | Rabbit  |
| Applications          | WB,IHC-p,IF,ELISA   |
| Species Reactivity    | Human,Mouse,Rat   |
| Recommended Dilutions |   |
| Immunogen             |   |
| Species               | Rabbit  |
| Storage               | -20°C/1 year  |
| Isotype               |   |
| Clonality             |   |
| Concentration         | 1 mg/ml   |
| Observed band         | 40kDa   |
| GenID?Human?          | ARFIP1  |
| Human Swiss-Prot No.  |   |
| Cellular localization |   |
| Alternative Names     | ARFIP1; Arfaptin-1; ADP-ribosylation factor-interacting protein 1   |
| Background            | function:Putative target protein of ADP-ribosylation factor.,similarity:Contains 1 AH domain.,subunit:Interacts with non-myristoylated GTP-bound ARF3, but not to GDP-bound ARF3, and also to ARF1. Binds with lower affinity to ARF5 and with very little affinity to ARF6.,tissue specificity:Ubiquitously expressed. Higher levels in liver, pancreas, placenta, skeletal muscle and heart., |