

GRP78/Bip Mouse mAb

Cat No: HR1AM2321

For research use only

Overview

Product Name	GRP78/Bip Mouse mAb
Source	Mouse
Applications	WB, IHC
Species Reactivity	Human,Rat
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:100-200
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG1 UEIDU
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	78kDakDa
GeneID?Human?	3309
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	HSPA5, GRP78,MIF2, BIP, Heat shock 70 kda protein 5
Background	Belongs to the heat shock protein 70 family. When protein folding is disturbed inside ER, BiP synthesis is increased.