




ApoL5 Polyclonal Antibody

Cat No: HR1AP8736

For research use only

Overview

Product Name	ApoL5 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GenID?Human?	APOL5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	APOL5; Apolipoprotein L5; Apolipoprotein L-V; ApoL-V
Background	apolipoprotein L5(APOL5) Homo sapiens This gene is a member of the apolipoprotein L gene family. The encoded protein is found in the cytoplasm, where it may affect the movement of lipids or allow the binding of lipids to organelles. [provided by RefSeq, Jul 2008],



bioelsa