

GPR158 Polyclonal Antibody

Cat No: HR1AP8023

For research use only

Overview

Product Name	GPR158 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	150kDa
GeneID?Human?	GPR158
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPR158; KIAA1136; Probable G-protein coupled receptor 158
Background	function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 3 family.,