

# EGFR (phospho Thr693) Polyclonal Antibody

**Cat No: HR1AP2411**

For research use only

## Overview

|                       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | EGFR (phospho Thr693) Polyclonal Antibody                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Source                | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Applications          | WB,IHC-p,IF,ELISA                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Species Reactivity    | Human,Mouse,Rat                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Recommended Dilutions |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Immunogen             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Species               | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Isotype               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Clonality             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Observed band         | 1,80,134kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| GeneID?Human?         | EGFR                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Cellular localization |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Alternative Names     | EGFR; ERBB; ERBB1; HER1; Epidermal growth factor receptor; Proto-oncogene c-ErbB-1; Receptor tyrosine-protein kinase erbB-1                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Background            | epidermal growth factor receptor(EGFR) Homo sapiens The protein encoded by this gene is a transmembrane glycoprotein that is a member of the protein kinase superfamily. This protein is a receptor for members of the epidermal growth factor family. EGFR is a cell surface protein that binds to epidermal growth factor. Binding of the protein to a ligand induces receptor dimerization and tyrosine autophosphorylation and leads to cell proliferation. Mutations in this gene are associated with lung cancer. [provided by RefSeq, Jun 2016], |