

GPR171 Polyclonal Antibody

Cat No: HR1AP6809

For research use only

Overview

Product Name	GPR171 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	GPR171
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GPR171; H963; Probable G-protein coupled receptor 171; G-protein coupled receptor H963
Background	function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 1 family.,