



BET3 Polyclonal Antibody

Cat No: HR1AP2879

For research use only

Overview

| | |
|-----------------------|---|
| Product Name | BET3 Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,IHC-p,ELISA |
| Species Reactivity | Human,Mouse,Rat |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 22kDa |
| GenID?Human? | TRAPPC3 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | TRAPPC3; BET3; CDABP0066; Trafficking protein particle complex subunit 3; BET3 homolog |
| Background | trafficking protein particle complex 3(TRAPPC3) <i>Homo sapiens</i> This gene encodes a component of the trafficking protein particle complex, which tethers transport vesicles to the cis-Golgi membrane. The encoded protein participates in the regulation of transport from the endoplasmic reticulum to the Golgi apparatus. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2012], |