

ApoOL Polyclonal Antibody

Cat No: HR1AP5927

For research use only

Overview

Product Name	ApoOL Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	30kDa
GeneID?Human?	APOOL
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	APOOL; CXorf33; FAM121A; Apolipoprotein O-like; Protein FAM121A
Background	apolipoprotein O like(APOOL) Homo sapiens This gene encodes a protein which contains an apolipoprotein O superfamily domain. This domain is found on proteins in circulating lipoprotein complexes. [provided by RefSeq, Sep 2011],